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
Established 1947

**Acquired Autoclave Control
and Engineering Ltd in 2001**

60 Employees in Group

Group Turnover of £4.5m

Nationwide Sales & Service Network

Based in Oldham, Lancashire

Manufacturer & Supplier of

Autoclaves & Thermal Products for Laboratory and Process Use

Through Strategic Partnerships

LTE Export to Every Continent

LTE Products include

Autoclaves


Ovens

Microbiological Incubators

Cooled Incubators

Freeze-Dryers

Environmental Rooms

General & Related Knowledge of Autoclave

Protozoa

Virus

Microorganism

Bacteria

Prion

Fungi

Bacteria

- ☐ Vegetative Form
- ☐ Reproductive Form (Spore)

Disinfection & Sterilization

Disinfection

By

Chemical

Disinfection

By

Heat

What is A_0 ?

$A_0 = \text{Second at } 80\text{ }^{\circ}\text{C}$

A_0	=	600	
70 °C	=	100	Min.
** 80 °C	=	10	Min.**
90 °C	=	1	Min.
93 °C	=	30	Sec.

Recommendation

☐ Human Waste Container ; $A_0 = 60$

☐ Surgical Instrument ; $A_0 = 600$

Mycobacteria ,Fungi and
Heat Sensitive Viruses } $A_0 = 600$

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
BHT INNOVA WD

BHT manufactures
small and simple as
well as large and
complex machines ...

E2

T5

Sterilization

Radiation

Chemical

Wet Heat

Gas

Dry Heat


Wet Heat

or


Moist Heat

Nature of Water


Solid



Heat

Liquid

Vapour

Heat

Steam

Pressure (lb / in²)

Temp (°C)

0

100

15

121

27



132

Steam

Contact

Target

Target

☐ Dirty

☐ Air

Sterile

Must Be

Clean

Air Removal

Autoclave Principle

- ☐ Sterile by Moisture Heat. (Steam)
- ☐ Increase Steam Temp. by Increase Steam Pressure.
- ☐ Require the time for Sterile.

Sterilizing Cycle

- ☐ Air Removal or Purging
- ☐ Sterilizing
- ☐ Cooling

General Knowledge of Autoclave

What is an Autoclave ?

A sterilising device using
high pressure steam


What sort of vessel ?

- Vertical cylindrical vessel
- Horizontal cylindrical vessel
- Horizontal rectangular vessel

Cylindrical Pressure Vessel

Vertical cylindrical vessel

□ Advantages

- Compact
- Easy to manufacture
- Small footprint
- Inexpensive

□ Disadvantages

- Difficult to use
- In chamber steam generation
- Normally only suitable for simple cycles
- Temperature not very stable

Vertical vessel autoclave

Horizontal cylindrical vessel

- Advantages
 - Ideal pressure vessel
 - Relatively inexpensive
 - Easy to manufacture
- Disadvantages
 - Doesn't maximise load capacity
 - Large footprint


Horizontal cylindrical vessel autoclave

Rectangular Pressure Vessel

Horizontal rectangular vessel

- Advantages
 - Maximises load capacity
 - Good steam and temperature distribution
 - Same footprint as cylindrical chambers
- Disadvantages
 - Requires thicker material
 - Not an ideal pressure container
 - More stringent welding requirements

Horizontal rectangular vessel autoclave

Materials ?

- Stainless steel
- What do we mean ??

Types of stainless steel

- 304
- 316 & 316 L
- 321

316L for GMP Autoclaves ?

- Very low carbon content
- Resistance to pitting
- Low tendency to rouging
- Reasonably easy to weld

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
Chemical Indicator ; MESA, Reven, USA

- Steam Class 4, ProChem SSI
- Steam Class 5, ProChem SSW
- Steam Class 5, ProChem EXT
- Bowie Dick Test, CI-BD110
(Vacuum Test)

Physical Indicator; Marathon, USA

- EZ Logger High Temp

- -40 °C to +125 °C

- SS

- 316 grade

- USB interface

- Compact size

- dia = 1.9 cm

Autoclave ; Certificate of Calibration (ISO/IEC 17025 : 2005)

105°C - 135°C ; BS 2646 Part 5 (1993)

$$\text{Survival} = 10^2$$

$$D_{121} \text{ Value} = 1.5$$

$$\text{Survival time} = 4.5 \text{ Minutes } (10^5)$$

$$\text{Survival time} = 6.0 \text{ Minutes } (10^6)$$

Cooling Cycle

- ☐ Fan
 - ☐ Water
 - ☐ Compressed Air
 - ☐ Vacuum

Reduce time

At least 50%

Drying

Vacuum (Decrease Pressure)

Hot Chamber

Type of Autoclave

Single Wall Autoclave ●

- Basic & General Autoclave
 - glassware, containers, instruments
- Liquid & Media Autoclave (with air ballast)
 - fluid in bottle/bag, media
- Food (F_0) & Low Temp ($< 121^{\circ}\text{C}$) Autoclave

LTE Scientific Ltd Touchclave-R

The Dynamic New Cylindrical Autoclave Range

Typical Loading Capacities

	TCR/40/H	TCR/60/H	TCR/60/V	TCR/120/ H	TCR/120/ V	TCR/160/ H	TCR/160/ V
1/2 litre bottle	14	21	18	48	42	63	63
1 litre bottle	9	14	14	23	32	30	48
2 litre bottle	5	8	6	13	14	18	14
1 litre flask	2	4	8	13	16	18	24

Why Will It Appeal To Customers?

- Applications
- Ease of Installation
- Ease of Operation
- Speed
- Price

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Thermal Data Printer (Optional)

Data Archiving (Optional)
5000 cycle

Temperature And Pressure Monitoring

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




Touchclave V + Dry

50 L, 80 L, 150 L

Complete with Touchscreen control panel
for Real Time of Pressure, Time and stage,
Temp. probe, Vacuum, Safety feature


Touchclave V + Dry

Popular Option

Data Printer – Details of program, temperature, pressure, time and stage.

Cat No. **AVTDP**


Touchclave V + Dry

Touchscreen Control Panel


1. Graphical representation of real time process
2. Start/Stop.
3. Running program.
4. Pressure indicator.
5. Temperature indicator.
6. Cycle parameters.

Data Archiving

- For long-term storage of sterilization process records
- Up to 5000 cycle profiles on a single flash card
- Can be printed off at any time


Touchscreen Control System

- Touchscreen
- 8-programs
- Pre-set and adjustable templates
- Alpha-numeric inputs
- Comprehensive alarms
- PIN security levels
- View real-time cycle or print last cycle


Door Closure (Push-n-Seal Door) and Pneumatic Door Seal

- Simply push the door closed
- Start the cycle
- Door locks automatically
- No levers, no bolts
- At the end of the cycle, the door automatically unlocks
- Unique in the market


Keep Warm Facility

- Control temperature of chamber at and of cycle at 65° c
- Allow the processed media in the chamber to remain liquid
- No need of extra water bath


Temperature and Pressure Monitoring

- Controlling by pressure provides greater temperature stability than controlling by temperature
- Important for certain media which are sensitive to temperature overcooked.
- Typical accuracy is $\pm 0.5^{\circ} \text{ c}$


Dynamic In-Chamber Cooling

- Two methods of cooling
- In-chamber cooling pushes air through the chamber via fine air inlet filter
- In chamber cooling can assist in drying certain load types
- Helps to speed up cooling times
- Unique feature


Load Sensed Process Timer

- Armour-flex probe
- Tougher than plastic or stripped-end probes
- Does not begin sterilization timer until load and drain reaches sterilizing temperature
- Does not release door until load is at safe temperature


Air Ballast

- Used for Media cycles
- Speeds up cooling process
- Maintains an over-pressure in chamber
- Compressed air enters chamber through a fine air filter
- No-one else offers this as a standard feature

Thermal Data Printer

LTE SCIENTIFIC

L7132

TCR

Report :- Cycle Setup
Program :- MEDIA 121:15 BALLAST
Date :- 18/09/07
Time :- 16:25

User Level For Start None
User Level For Approve None
Door Safety Temp 000.0 DegC
Print During Cycle YES
Discard Cycle NO
Batch Code Required NO
Cooling O/Ride Fails Cycle NO
Keep Warm YES

DOOR SEALING
Monitor Time 0003 Mins

STEAM PURGE
Monitor Time 0120 Mins
Set Temp 090.0 DegC
Set Time 0600 Secs
Safety Pressure 2000 mBar

HEATUP
Monitor Time 0120 Mins
Boost Factor None
No Boost Band 000.0 DegC

STABILISATION

Monitor Time 0010 Mins
Stabilisation Time 0060 Secs

STERILIZATION

Monitor Time 0060 Mins
Sterilize Time 0900 Secs
Sterilize Pressure 0000 mBar
Sterilize Temp 121.0 DegC
Dropout Temp Limit 001.0 DegC
Overshoot Temp Limit 005.0 DegC
Logging Interval 0060 Secs

COOLING

Cooling Type Ballast
Cooling Temp 080.0 DegC
Minimum Cooling Time 0000 Mins
Cooling Assist Delay Pressure 0000 mBar
Ballast Low Pressure 1900 mBar
Ballast High Pressure 2000 mBar
Final Discharge Pressure 1150 mBar
Final Discharge Path Exhaust
Cooling Interlock NO
Logging Interval 0300 Secs

AIR ADMISSION

Monitor Time 0010 Mins
Settling Time 0030 Secs

LTE SCIENTIFIC
GREENBRIDGE LANE
GREENFIELD
OLDHAM
LTE SERVICE CALL
Tel:01457-876221

LTE SCIENTIFIC

L7132

TCR

Report :- Cycle
Program :- MEDIA 121:15 DYNAMIC
Operator :- N/A
Cycle No. :- 000023
Date :- 19/09/07
Start Time:- 09:27
Steri Temp:- 121.0 oC
Steri Time:- 00:15:00

Time	Cham	Vent	Load
H:M:S	mBar	oC	oC

DOOR SEALING

00:00:00	0986	027.0	025.7
00:00:08	0952	026.9	025.7

STEAM PURGE

00:00:08	0952	026.9	025.7
00:20:16	1452	110.4	106.2

HEATUP

00:20:16	1452	110.4	106.2
00:28:57	2085	121.8	121.0

STABILISATION

00:28:57	2085	121.8	121.0
00:29:57	2094	122.1	121.2

STERILIZATION

00:29:57	2094	122.1	121.2
00:30:57	2075	121.8	121.3
00:31:57	2102	122.1	121.3
00:32:57	2077	121.9	121.4
00:33:57	2107	122.0	121.4
00:34:57	2092	122.1	121.5
00:35:57	2085	121.8	121.5
00:36:57	2108	122.2	121.5
00:37:57	2089	122.0	121.5
00:38:57	2087	121.8	121.6
00:39:57	2114	122.1	121.6
00:40:57	2106	122.1	121.7
00:41:57	2093	122.0	121.7
00:42:57	2080	121.8	121.7
00:43:57	2111	121.9	121.7
00:44:57	2113	122.1	121.8
00:44:57	2113	122.1	121.8

COOLING

00:44:57	2113	122.1	121.8
01:19:41	1080	057.2	080.0

AIR ADMISSION

01:19:41	1080	057.2	080.0
01:20:14	1015	052.9	079.9

F0 Achieved :- 00:22:24

Cycle No. :- 000023
Errors :- No
Warnings :- No
Result :- Passed
Cycle Time:- 01:20:14
Approved :- N/A

Cycle Complete

Double Wall Autoclave ●

Applications

- Central Sterile Services Department (CSSD)
- Production → ● Laboratory

Features

- Capacities 300 to 10,000 litres (10 – 350 cu.ft)
- Wide range of applications
- Pressure vessels compliant with and independently tested to meet European and international safety standards
- All chambers are manufactured from high *grade 316L stainless steel*

Pressure Equipment Regulations 1999 *Certificate of Conformity- EC Product Verification* Form PER F

CERTIFICATE REFERENCE: CEN-005517/001

MANUFACTURER NAME	LTE Scientific Ltd
ADDRESS	Greenbridge Lane, Greenfield, Oldham OL3 7EN
MANUFACTURER REFERENCE	J040101 (Mr. Simon Parr)

EQUIPMENT DESCRIPTION	Autoclave SCF36MP
EC TYPE/DESIGN	SS42071113
EXAMINATION REFERENCE	
CONSTRUCTION CODE	PD5500:2000 CAT 3
SERIAL NUMBER	J040101
DRAWING NUMBER(S)	50116 rev 4 & 50115 rev 2
DESIGN PARAMETERS	Chamber 2.6 barg / FV @150 °C Jacket 2.6 barg @ 150 °C
PRESSURE TEST DETAILS	Chamber 5.2 bar water 30 min on 09/03/2004 Jacket 6.2 bar water 30 min on 09/03/2004

This is to certify that the pressure equipment identified above was subjected to EC Product Verification, during and upon completion of manufacture, as required by the code of construction and inspection specification, and was found to meet with the provisions of the Pressure Equipment Regulations, subject to any conditions of use specified below.

Zurich Risk Services
Authorised Signatory

Title: Senior Engineer

Notified Body No 0037
Date: 18 March 2004

Lei Zhu

PC ZCERT-4 10/03

Form X (concluded)

Vessel Serial

9904

Certificate of Compliance (design)

We hereby certify that the design of this vessel complies with PD 5500 :2000

For manufacturer (see note) 6):

Date:

Position: **DESIGN ENGINEER.**

Name of company: **LTE SCIENTIFIC LTD**

We hereby confirm that we have checked the design of the above vessel and that this complies with PD 5500:2000

For Inspecting Authority:

Date:

Position: **SENIOR ENGINEER**

Name of company: **ZURICH.**

Certificate of Compliance (construction and testing)

Vessel Serial

9904

We hereby certify that this vessel has been constructed and tested in compliance with PD5500:2000

For manufacturer:

Date:

10/03/2004

Position: **DESIGN ENGINEER**

Name of company: **M H VESSELS & FABRICATIONS LTD**

We hereby confirm that the construction and testing of the vessel has been carried out under our surveillance and that to the best of our knowledge and belief all aspects of this work comply with PD 5500

For Inspecting Authority:

Date:

10/03/2004

Position: **ENGINEER SURVEYOR** Name of company:

ZURICH

NOTE 1. The suffix 'XX' is to be added to the serial number of each vessel for which any deviations or concessions have been authorised (see 1.5.2.2 (e)).

NOTE 2. The design conditions associated with the operational duties specified by the purchaser should be give

If a purchaser wishes to change the operational duty of a vessel, revised design conditions, consistent with the vessel scantlings, will be established separately, as appropriate.

NOTE 3. Where the covers operation below 0 C the various combinations of temperature ,pressure and calculated membrane stresses considered in determining the design minimum temperature (see 3.2.5 and annex D) should be stated.

NOTE 4. Where appropriate, cross reference to drawing or specifications will suffice.

NOTE 5. Where a vessel is tested in a different orientation to that in which it will normally operate, this should be stated.

NOTE 6. This part of the Certificate to be signed by the purchaser in cases where the purchaser elects to perform the design function (see 1.4.4).

Features (Cont.)

- **Emergency stop button for manual shutdown**
- **LTE sterilizers are equipped with our advanced touchclave touchscreen control system**
- **Internal data archiving ; 5,000 cycles**
- **Integral thermal or ribbon printer included**

Combined vacuum and air ballast

Sterling Anti-Bacterial powder coating

The Door

- **316L SS**
- **Pneumatically actuated for smooth movement**
- **Interlocked system in a double door model**
- **Horizontal or vertical sliding door**

The Door

- “Easy-glide” movement.
- Floating door seal, using built-in compressor to push seal onto door and then retract at end of cycle.
- High level of reliability
- Easy-to-replace

Touchclave Control System

- State-of-the-art Omron PLC
- IP rated Omron NT31 Touchscreen

“Cycle Select” Screen


Touchclave Control System

Delayed start function

Touchclave Control System

Page 1 of alarm screen

Data Archiving

- 5000 cycles could be stored on a single flash card
 - Option : download data to a PC

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

0.00


Process in Autoclave

- **Air Removal phase**
- **Sterilization phase**
- **Cooling / Drying phase**


Conclusion

- Feature of Jacket Type
- (Double wall)
- Feature of vacuum

Result

- Load comes out as
- Completely dry & Complete sterilization
- (discard infectious garbage)

Technical and Ordering Information

Model	Capacity Litres / Cu.Ft.	Single Entry (S/E) Pass Through (P/T)	Door Closure	Chamber H x W x D (mm)	Overall Dimensions H x W x D (mm)	Approx. Weight (Kg)
SCF 10	326 / 11	S/E	HORIZONTAL	660 x 660 x 750	1762 x 1880 x 1690	800
UCF10	326 / 11	S/E	VERTICAL		1800 x 1290 x 1626	800
2SCF10	326 / 11	P/T	HORIZONTAL		1762 x 1880 x 1280	800
2UCF10	326 / 11	P/T	VERTICAL		1800 x 1290 x 1152	800
SCF15	435 / 15	S/E	HORIZONTAL	660 x 660 x 1000	1762 x 1880 x 1940	875
UCF15	435 / 15	S/E	VERTICAL		1800 x 1290 x 1850	875
2SCF15	435 / 15	P/T	HORIZONTAL		1762 x 1880 x 1530	875
2UCF15	435 / 15	P/T	VERTICAL		1800 x 1290 x 1402	875
UCF20	566 / 20	S/E	VERTICAL	660 x 660 x 1250	1800 x 1250 x 2130	950
SCF21	636 / 21	S/E	HORIZONTAL	965 x 660 x 1000	1762 x 1880 x 2190	950
2SCF21	636 / 21	P/T	HORIZONTAL		1762 x 1880 x 1780	950
SCF28	796 / 28	S/E	HORIZONTAL	965 x 660 x 1250	1762 x 1880 x 2390	1250
2SCF28	796 / 28	P/T	HORIZONTAL		1762 x 1880 x 1980	1250
SCF36	1052 / 36	S/E	HORIZONTAL	1100 x 660 x 1450	1762 x 1880 x 2390	1500
2SCF36	1052 / 36	P/T	HORIZONTAL		1762 x 1880 x 1980	1500

Note: In addition to the product dimensions listed above, LTE Scientific is also able to supply autoclaves to your precise requirements.

	10/11 cu.ft.	15 cu.ft.	20/21 cu.ft.	28 cu.ft.	36 cu.ft.
Steam Flow Rate – Average (kg/hr)	26	30	40	45	50
Steam Flow Rate – Maximum (kg/hr)	52	65	80	90	100
Steam Consumption per Cycle (kg – 30 min cycle)	13	15	20	23	25
Condensate Flow Rate – Maximum (kg/hr)	60				
Condensate Discharge/Drain Requirements	One x 100mm (4”) diameter floor gully capable of withstanding temperatures in excess of 110°C				
Steam Supply (where applicable)	One x 25mm (1”) diameter steam supply to 5 bar				
Water Flow Rate – Maximum (l/min)	7.6				
Water Supply Pressure (bar)	2bar (flow)				
Drain Size (mm)	54				
Drain Type	Sealed				
Drain Internal Type	Vented				
Vent Size and Type	35mm Sealed				
Safety Valve Outlet Size (mm)	35				
Compressed Air Flow Rate (m ³)	0.5				
Compressed Air Supply Pressure (bar)	4-8				
Compressed Air Consumption Per Cycle (m ³)	1.0				
Electricity Voltage	415V 3-Phase Neutral + Earth				
Electricity Current (Excl. Steam Generator)	20 amps per phase				
Electricity Current For Steam Generator (if applicable) Max. (amps per phase)	65	80	80-90	80-90	100-130
Electricity Max. Power – Direct Steam Models (kw)	5				
Electricity Max. Power – Steam Generators (kw)	40	45	45-55	45-55	60-70
Air Filter (Air Removal) Expected Life	6 months				
Sound Power	As specified in EN ISO 3746:1996				

Alarm list

13. ALARMS

13.1 The Alarms System

All alarms remain triggered until reset either by the reset button on the alarm screen or by completing the cycle.

It is impossible to clear alarms in the middle of a cycle.

An alarm will fail to clear if the condition causing the alarm remains.

Below is a list of alarms and their descriptions:

13.2 Alarm List

Error	Description
Cycle Aborted	A User has manually aborted the cycle from the option in the main menu
Emergency Stop Activated	The emergency stop button has been activated and the reset button has not been reset.
Steam Failure	No Steam to the machine.
Water Failure	No Water to the machine.
Air Failure	No Air to the machine.
Chamber Transducer Failure	Chamber Transducer out-of-range or not connected.
Jacket Transducer Failure	Jacket Transducer out-of-range or not connected.
Drain Probe Failure	Drain Probe out-of-range or not connected.
Air Detector Probe Failure	Air Detector Probe out-of-range or not connected.
Vent Probe Failure	Vent Probe out-of-range or not connected.
Load Probe Failure	Load Probe out-of-range or not connected.
Filter Probe Failure	Filter Probe out-of-range or not connected.
Extra Probe Failure	Extra Probe out-of-range or not connected.
Door 1 Seal Failure	No signal is being received from the seal/clamp sensors for door 1 when the door should be sealed/clamped.
Door 2 Seal Failure	No signal is being received from the seal/clamp sensors for door 2 when the door should be sealed/clamped.
Power Failure	The machine has been turned off mid-cycle
Archive Save Failure	An out of cycle error that occurs when the last cycle has failed to save to the flash card. It will however try every 30 seconds until successful.
Unsafe Probe Temps	The Door Safety Temp has been exceeded between a cycle finishing and a door being opened.
Door 1 Not Closed	No signal is being received from the closed sensor for door 1.
Door 2 Not Closed	No signal is being received from the closed sensor for door 2.
Door 1 Not Bolted	No signal is being received from the bolted sensor for door 1.
Door 2 Not Bolted	No signal is being received from the bolted sensor for door 2.
Steam Gen Timeout	Stage time has exceeded the stage monitor time in the Awaiting Steam Gen stage.
Door Sealing Timeout	Stage time has exceeded the stage monitor time in the Door Sealing stage.
Pre Vacuum Timeout	Stage time has exceeded the stage monitor time in the Pre Vacuum stage.
Free Steaming Timeout	Stage time has exceeded the stage monitor time in the Free Steaming stage.

Operator Manual

Steam Flushing Timeout	Stage time has exceeded the stage monitor time in the Steam Flushing stage.
Negative Pulse Timeout	Stage time has exceeded the stage monitor time in the Negative Pulsing stage.
Air Detector Test Timeout	Stage time has exceeded the stage monitor time in the Air Detector Test stage.
Positive Pulse Timeout	Stage time has exceeded the stage monitor time in the Positive Pulsing stage.
Heating Timeout	Stage time has exceeded the stage monitor time in the Heating stage.
Stabilisation Timeout	Stage time has exceeded the stage monitor time in the Stabilisation stage.
Sterilisation Timeout	Stage time has exceeded the stage monitor time in the Sterilisation stage.
Leak Vacuum Timeout	Stage time has exceeded the stage monitor time in the Leak Vacuum stage.
Leak Stabilisation Timeout	Stage time has exceeded the stage monitor time in the Leak Stabilisation stage.
Leak Test Timeout	Stage time has exceeded the stage monitor time in the Leak Test stage.
Air Detector Fail	The air detector test detected air in process cycle during the Air Detector Test stage.
Leak out of Test Range	The pressure has risen too much in stabilisation phase and the test cannot be performed because it is above a set value.
Sterilise Under Temp	A Control Probe has dropped more than the Sterilise Dropout Temperature lower than the Sterilise Temperature.
Sterilise Over Temp	A Control Probe has gone more than the Sterilise Over Temperature above than the Sterilise Temperature.
Sterilise Under Pressure	The Chamber Pressure has dropped more than the Sterilise Under Pressure lower than the Sterilise Pressure.
Sterilise Over Pressure	The Chamber Pressure has gone more than the Sterilise Over Pressure above than the Sterilise Pressure.
Discard Hold Timeout	Stage time has exceeded the stage monitor time in the Discard Hold stage.
Cooling Timeout	Stage time has exceeded the stage monitor time in the Cooling stage.
Draining Timeout	Stage time has exceeded the stage monitor time in the Draining stage.
Drying Timeout	Stage time has exceeded the stage monitor time in the Drying stage.
Airwash Timeout	Stage time has exceeded the stage monitor time in the Airwash stage.
Air Admission Timeout	Stage time has exceeded the stage monitor time in the Air Admission stage.
Jacket Wait Timeout	Stage time has exceeded the stage monitor time in the Awaiting Jacket stage.
Water Fill Timeout	Stage time has exceeded the stage monitor time in the Water Fill stage.
Clean Steam Failure	No Clean Steam to the machine.
Cooling Override Used	Depending on whether using cooling override is set up to fault the cycle this error indicates that cooling override was used.
Steam Gen Water Low	Steam Gen Water at Extremely Low Levels
Steam Gen Over Temp	Heater Element Cutout
Engineering Override Used	A user has used an option in the engineer override mode, which could affect the integrity of the cycle.

LTE Scientific

Autoclaves and Thermal Products

In the past 25 years, L T E Scientific have supplied and installed over 1000 Autoclaves and associated products. An extract from our customer reference list is as follows :

Chinese University of Hong Kong CUHK – Hong Kong

National Institutes of Health, Bangkok, Thailand

Central Fisheries Unit, ICAR, India

Chulalongkorn University, Bangkok, Thailand

Hong Kong University HKU – Hong Kong

Universiti Putra Malaya, UPM, Malaysia

Islandwide Scientific, Colombo, Sri Lanka

Rizal Pharmaceutical, Kuala Lumpur, Malaysia

Development Centre for Biotechnology, Taipei, Taiwan

Technimex, Hanoi, Vietnam

Johnson & Johnson, Malaysia

National University of Singapore, Singapore

Sanofi Winthrop, UK

Public Health Laboratory, Kowloon, Hong Kong

Glaxo Smith Kline, UK

South East Asia Medical Epidemiology Organisation, Bangkok, Thailand

Fortune Electronics, China

Rolls Royce, Derby, UK

Merck Sharp & Dohme, UK

Defence Services Organisation, Singapore

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
Advantage Point & Benefit of Touchclave Lab 360 & 450 KE

- Large Chamber Size
- Dry (hot air oven is unnecessary)
 - comfortable & decrease contaminated risk
 - Saving Energy & Green Policy
- Complete sterilization

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
CO₂ EMISSION

**From Waste
Incinerator**

CO₂ EMISSION From Waste Incinerator

- Waste : MSW = Municipal Solid Waste
- Emission depend on composition of the waste (Carbon Compound)
- Emission Data is from the associated website

Typical CO₂ emission is 0.7-1.7 ton /ton of waste

From : <https://zerowasteeurope.eu/2020/03/understanding-the-carbon-impacts-of-waste-to-energy/>

The CO₂ emitted per tonne of waste incinerated depends on the composition of the waste. Eurostat statistics show that the majority of the increase in waste incinerated comes from so-called residual municipal solid waste (MSW) which produces significant amounts of CO₂.

Each tonne of MSW incinerated typically releases between 0.7 and 1.7 tonnes of CO₂ [3]. This includes emissions of both fossil CO₂ (e.g. from burning plastics) and biogenic CO₂ (e.g. from burning wood, paper and food). Although biogenic CO₂ is directly released into the atmosphere making a significant contribution to climate change, only the CO₂ emissions from fossil sources will be considered for the purposes of a global analysis – an important loophole in GHG emissions accountability.

Typical CO₂ emission is 0.7-1.2 ton /ton of waste

From : https://www.ipcc-nggip.iges.or.jp/public/gp/bgp/5_3_Waste_Incineration.pdf

Waste_Incineration.pdf

๑๖

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ABSTRACT

The incineration of municipal waste involves the generation of climate-relevant emissions. These are mainly emissions of CO₂ (carbon dioxide) as well as N₂O (nitrous oxide), NO_x (oxides of nitrogen) NH₃ (ammonia) and organic C, measured as total carbon. CH₄ (methane) is not generated in waste incineration during normal operation. It only arises in particular, exceptional, cases and to a small extent (from waste remaining in the waste bunker), so that in quantitative terms CH₄ is not to be regarded as climate-relevant. CO₂ constitutes the chief climate-relevant emission of waste incineration and is considerably higher, by not less than 10², than the other emissions.

Formulas (1) and (2) are to be used for the purpose of compiling an inventory of greenhouse gas emissions, taking the following into account:

The incineration of 1 Mg of municipal waste in MSW incinerators is associated with the production/release of about 0.7 to 1.2 Mg of carbon dioxide (CO₂ output). The proportion of carbon of biogenic origin is usually in the range of 33 to 50 percent. The climate-relevant CO₂ emissions from waste incineration are determined by the proportion of waste whose carbon compounds are assumed to be of fossil origin. The allocation to fossil or biogenic carbon has a crucial influence on the calculated amounts of climate-relevant CO₂ emissions.

Annex 1 contains a description of a method to calculate the energy credit for the use of waste as a substitute for fossil fuel (MSW incineration plants with energy recovery). Formulas (3) and (4) are to be used for the purpose of a comparative evaluation of the climate-relevant emissions from waste incineration in relation to those of other types of energy production. A factor that has a decisive influence on the calculated amounts of climate-relevant emissions from waste incineration plants with energy utilisation is the credit allowed or allowable due to the substitution of energy from fossil fuels. The latter in turn is influenced by the energy carriers used as a basis to calculate the emission factor of the power plant mix.

CO₂ EMISSION From Waste Incinerator

- Approx. 1 ton per ton of MSW
- **Zero Waste Europe** is the European network of communities, local leaders, experts, and change agents working towards the elimination of waste in our society.
- **The Intergovernmental Panel on Climate Change (IPCC)** is the United Nations body for assessing the science related to climate change.

Thank you